













Description		Order Code
Cleaning Holder		PD-CTH
Cleaning Adaptor Luer M		PD-CTAD-LM
Cleaning Adaptor Luer F		PD-CTAD-LF
Cleaning Adaptor SM		PD-CTAD-SM
Cleaning Adaptor 1T		PD-CTAD-1T
Cleaning Adaptor IN		PD-CTAD-IN
Cleaning Adaptor HM		PD-CTAD-HM
Cleaning Adaptor WPD		PD-CTAD-WPD
Cleaning Adaptor G		PD-CTAD-G
Cleaning Tube 6.0/L650		PD-CTTU-606501
Cleaning Tube 6.0/L50		PD-CTTU-600504
Cleaning Tube 1.5/L25		PD-CTTU-150252
EZ Spray Adaptor		PD-EZ-M
Spray Adaptor for Slim		P200-EZ-SM
Spray Adaptor for Bone Saw		P200-EZ-BM
Spray Adaptor for WPD		P200-EZ-WPD
PANA SPRAY Plus		Z182100

Bundled Item: PD-CTTU-606501
PD-CTTU-600504

Bundled Item: PD-CTTU-150252

Cleaning Tool

for Primado2 Attachments



Primado2

Efficient cleaning of contaminants in the lumen

Reduce the risk of excessive heat development and malfunction of attachments due to insufficient cleaning




Contaminants in the lumen

Specially designed Cleaning Holder


Maintains balance of flow rate in each outlet

Comparison table of flow rate*

- Controlled Design
- Uncontrolled Design



*Based on the result of in-house testing



Choice of cleaning options



Manual cleaning



Washer disinfector



Cleaning process steps for manual cleaning



*Please be sure to read the reprocessing manual.
 *These Cleaning Tools can not be used with the Bone Saw Attachment.
 *Cleaning Adaptor IN can not be used for automated cleaning.